


Add: No.2000 Shunhua Rd, High-Tech Zone , Jinan City, Shandong Province, China Contact: +86-18653132120

Product name:Potassium Amylxanthate

CAS:2720-73-2

Molecular Weight:204.39



٠					٠
L	,		L	-	ı
	×.	. 1	r	-	t
		٤.			

Item	Specification		
Boiling point	497.18°C[at 101 325 Pa]		
Density	1.24[at 20°C]		
Vapor pressure	0Pa at 25°C		
Purity	97.0%		
Storage conditions	Inert atmosphere,Room Temperature		

Application:

Potassium Amylsantate is a strong collector mainly used in the flotation of non-ferrous metal minerals that require strong collection power without selectivity. For example, it is a good collector for flotation oxidized sulfide ore or oxidized copper ore and oxidized lead ore (sulfided with sodium sulfide or sodium hydrosulfide). This product can also achieve good separation effects on copper nickel sulfide ores and gold bearing pyrite flotation.

Package:

Usually packed in 25kg/drum, and also can be do customized package.